### **Supplemental Methods**

### Human adipose tissue collection

17 patients with symptomatic, persistent or long-standing persistent AF undergoing videoassisted thoracoscopic pulmonary vein isolation were recruited at the Amsterdam UMC, University of Amsterdam, between 2019 and 2022. Written informed consent was obtained from all participants<sup>1</sup>. The study was conducted in accordance with the Declaration of Helsinki as revised in 2013 and with the Medical Research Involving Human Subjects Act (WMO) and was approved by the Institutional Review Board of the Amsterdam AMC. AF was documented on ECG, Holter or pacemaker electrogram at least once in the 12 months preceding surgery. Patients with symptoms of severe heart failure or a left ventricular ejection fraction <35% were excluded from the study. The clinical characteristics of patients are listed in **Supplemental Table 1**. Patient EAT and SAT samples or the secretome derived from it were used for flow cytometry, miRNA sequencing and/or miRNA validation experiments (**Supplemental Table 1**).

During thoracoscopic ablation the left atrial appendage (LAA) was resected as standard procedure. The atrial tissue was collected and the EAT overlying the LAA was removed. A biopsy of approximately 1 cm<sup>3</sup> of SAT was obtained from one of the thoracoscopic entry sites. Adipose tissue was transported in phosphate-buffered saline to the laboratory for secretome extraction, or immediately snap-frozen for subsequent RNA isolation.

### Adipose tissue secretome collection

The method for secretome collection has been described previously<sup>5</sup>. Briefly, EAT and SAT were rinsed and minced into 1-2 mm<sup>3</sup> explants. Each piece was transferred to a well of a 96-

well plate and incubated in 100  $\mu$ L culture medium for 24 hours at 37°C and 5% CO<sub>2</sub> in order to stimulate adipokine and secretome release<sup>2</sup>. EAT and SAT secretome samples from each patient were collected, pooled, filtered (0.45  $\mu$ m pore size hydrophilic PVDF membrane, Millipore, SLHV033RB) and stored at -80°C. Adipose tissue was weighted after secretome collection to ensure similar adipose tissue weight per mL of secretome (**Supplemental Figure 1**).

#### **Concentration measurements of extracellular vesicles by flow cytometry**

Flow cytometry (A60-Micro, Apogee Flow Systems, Hemel Hempstead, UK) was used to determine the EV concentration<sup>3</sup> in adipose tissue secretome samples collected from 7 patients (clinical characteristics provided in **Supplemental Table 1**) and control culture medium without serum. For each of the 7 patients, we collected EAT and SAT samples, and extracted its secretome.

To achieve a count rate below 5000 events/seconds upon measuring, thereby preventing swarm detection<sup>4</sup>, adipose tissue conditioned culture medium samples were diluted 1.5 to 2-fold in Dulbecco's phosphate-buffered saline prior to antibody labelling. To label EVs, 20  $\mu$ L of diluted sample was incubated with 2.5  $\mu$ L of antibodies (antibody references are provided in MIFlowCyt-EV supplement) and kept in the dark at room temperature for 2 hours. We labelled EVs using tetraspanin markers (CD63 and CD9), and lactadherin to identify phosphatidylserine exposing EVs. To decrease background fluorescence by unbound reagents, 200  $\mu$ L of Dulbecco's phosphate-buffered saline was added to the sample prior to measurement. Samples were measured for 120 s at a flow rate of 3.01  $\mu$ L/min. The trigger threshold was set to 24 arbitrary units of the side scattering detector (side scattering cross section of 8 nm<sup>2</sup>). The reported concentration of EVs per mL of secretome refers to the number of particles (i) exceeding the side scattering threshold, (ii) having a diameter <1,000 nm, and (iii) exceeding the fluorescent gate corresponding to the used labels.

To improve the reproducibility, flow cytometry data were reported according to the MIFlowCyt-EV framework<sup>5</sup>. We calibrated all detectors, applied assay controls to confirm that signals originate from EVs, and automated data calibration and processing with custom-built software<sup>6</sup> (MATLAB R2020b, Mathworks, USA).

### Neonatal rat ventricular myocytes isolation

All animal experiments were approved by the local Animal Experiments Committee (Academic Medical Center, University of Amsterdam) and carried out in compliance with the Guide for the Care and Use of Laboratory Animals and in accordance with national and institutional guidelines. NRVMs were isolated from 1-to-2-day-old Wistar rats (Janvier labs) as previously described. Pups were anesthetized with isoflurane and hearts were excised after decapitation. Ventricular myocardium was cut into pieces and dissociated with trypsin (1 mg/mL;Sigma), and collagenase type 2 (1 mg/mL;Worthington). Cells were collected and resuspended in TUNG culture medium [M199 medium (Gibco) supplemented with 10% heat inactivated fetal bovine serum (FBS, Gibco), 1% HEPES (Gibco), 5000 U/L penicillin-G (Sigma), 2 mg/L vitamin B12 (Sigma-Aldrich), 3.5 g/L glucose, 1% non-essential amino acids (Gibco), 1% L-glutamine (Gibco)]. The cell suspension was pre-plated for 2 hours to minimize fibroblast contamination. The remaining myocardial cells were plated on microelectrode arrays (MEAs, Multi-Channel Systems MCS GmbH, Reutlingen, Germany) or 6-well-plates coated with fibronectin (Corning, Ref 356008, Bedford, MA, USA) at a density of 1.4x10<sup>5</sup> cells per cm<sup>2</sup>. NRVMs were cultured at 37°C and 5% CO<sub>2</sub> for 3 days in TUNG medium.

### Transfection with miRNA mimics

At 20–24 hours after seeding, NRVMs were transfected with 100 nM negative control siRNA (Qiagen, Ref 1027417, Germantown, MD, USA), 100 nM Syn-hsa-miR-133a-3p miScript

miRNA Mimic or 100 nM Syn-hsa-miR-1-3p miScript miRNA Mimic (Qiagen, Ref 219600, Germantown, MD, USA) using Lipofectamine<sup>™</sup> RNAiMAX Transfection Reagent (Thermo Fisher Scientific, Ref 13778150, Carlsbad, CA, USA) in M199 medium (Gibco) without supplements. After 6 hours of incubation, cell culture medium was replaced by TUNG with 2% FBS. Cells were harvested or measured at 48 hours after transfection. miRNA mimic sequences are found in **Supplemental Table 3**.

## qRT-PCR

RNA was extracted from EAT, SAT and NRVMs with TRI reagent (Sigma-Aldrich) following manufacturer's instructions.

Total RNA was converted into cDNA using the qScript microRNA cDNA Synthesis kit (Quantabio, CA, USA) according to the manufacturers' instructions.

Briefly, miRNAs were polyadenylated in a poly(A) polymerase reaction. qScript Reverse Transcriptase and other necessary reagents were added to convert the poly(A) tailed miRNAs and poly(A) containing mRNA into cDNA using an oligo-dT primer with an adapter. This miRNA-based cDNA was amplified in a RT-qPCR reaction with a miRNA-specific forward primer (sequences reported in **Supplemental Table 4**) and a general reverse primer binding to the adapter of the oligo-dT primer. RT-qPCR of mRNA was performed with gene-specific intron-spanning forward and reverse primers. RT-qPCR was then performed on a LightCycler 480 (Roche) using the PerfeCTa SYBR Green Kit (Quantabio, USA) and the following program:.

Results were analyzed using LinRegPCR software<sup>7</sup>. For *Kcnj2* gene expression quantification, the geometric mean of the reference genes *Hprt1* and *Eef1e1* was used to normalize expression levels and depicted gene expression levels are shown relative to the

control condition in each NRVM isolation. *U6* was used as a reference gene to normalize miRNA expression levels.

### Electrophysiology

Electrical mapping was performed as previously described. Briefly, 10 minutes before measurements, culture media was replaced by a modified Tyrode's solution (36.5°C) containing (mM): NaCl 140, KCl 5.4, CaCl<sub>2</sub> 1.8, MgCl<sub>2</sub> 1.0, glucose 5.5, HEPES 5.0; pH 7.4 (adjusted with NaOH).

NRVMs cultured on multi-electrode arrays (MEAs, 60 electrodes terminals, 60EcoMEA-Glass-gr, MultiChannel Systems MCS GmbH) were stimulated using a bipolar extracellular stimulus electrode (1.5x diastolic stimulation threshold, 2ms pulse). A reference electrode was placed inside the MEAs, not in contact with the cells. Unipolar electrograms were recorded simultaneously with a 256-channel amplifier (BioSemi, 24 bit dynamic range, 122.07 nV LSB, total noise 0.5  $\mu$ V). Signals were recorded with a sampling frequency of 2048 Hz (bandwidth (-3dB) DC–400 Hz).

Data analysis was performed using a custom-made data analysis program written in Matlab 2006b (The MathWorks Inc, Natick, MA<sup>8</sup>). Local activation time was determined as the interval from the stimulation artefact to the minimum derivative of the local QRS, and used to construct activation maps.

Blinded analysis was performed to measure conduction velocity, which was determined along lines perpendicular to isochronal lines by dividing the difference in local activation times by the electrode distance. Three conduction velocity measurements were performed and averaged per monolayer. Heterogeneity in conduction was measured based on the method described by Lammers et al<sup>9</sup>. Heterogeneity in conduction was determined as follows: at every electrode recording location, the maximal difference with neighboring activation times was calculated, known as the phase difference. These local phase differences were then visualized on a phase map, revealing the spatial pattern of heterogeneity in conduction. Subsequently, the median value of heterogeneity was derived from each of these maps.

### EAT secretome-derived EV isolation

Frozen EAT secretome samples from 8 patients were thawed on ice, and concentrated to 1 mL using Amicon Ultra-15mL Centrifugal Filter Units. The clinical characteristics of patients are listed in **Supplemental Table 1**.

EAT secretome-derived extracellular vesicles (EATS-EVs) were isolated using qEVoriginal 70 nm columns (IZON Science LTD) combined with an automatic fraction collector. The EV enriched fractions 2,3 and 4 (between 0.5 and 2 mL of the Purified Collection Volume) were pooled and concentrated to 200  $\mu$ L using Amicon Ultra-2mL centrifugal Filter Unit with a 10 kDa cut off (Merck#UFC201024). RNA from EATS-EVs was isolated using the miRNeasy serum/plasma kit (QIAgen) according to the manufacturers' protocol.

## Small RNA library preparation and sequencing

miRNA sequencing was performed by ExBiome as previously described<sup>10</sup> and as detailed below. Small RNA libraries were prepared using a previously published IsoSeek protocol<sup>10</sup>. Briefly, custom designed adapters were used in combination with the NEBNext Multiplex Small RNA Library Prep Kit for Illumina. The number of PCR amplification cycles was 20. Libraries were quantified with the KAPA Library Quantification Kit (Roche), mixed at equimolar concentration (1.5 nM each) and sequenced on a MiSeq platform (PE50, v2 300 cycles, with the addition of 10% PhiX).

### Processing of sequencing data and miRNA profiling

Pre-processing and mapping of adapter trimmed reads were performed using the latest version of sRNAbench command line tool<sup>11</sup>. Default parameters were used for all analysis steps after pre-processing and miRBase v22.1<sup>12</sup> was used as miRNA reference. Quality control of samples was carried out using mirnaQC<sup>13</sup> to rule out technical differences between libraries.

### **Prediction analysis**

GeneTrail2<sup>14</sup> was used to perform Gene Ontology (GO) analysis on 'cellular component' and 'biological process' terms based on predicted targets. To perform this analysis, the 156 miRNAs identified in all 8 patients were used as input in GeneTrail2. The indicated q-value corresponds to the adjusted *P*-value for multiple testing using the Benjamin-Yekutieli method. The indicated enrichment value corresponds to the number of hits divided by the expected score. TargetScan<sup>15</sup> was used to predict miRNAs targeting *KCNJ2*, using as input list the 100 most highly expressed miRNAs identified in all 8 patients. In addition, the miRDB<sup>16</sup> prediction database was used to identify potassium channel targets of miRNAs identified as regulators of resting membrane potential, but not predicted to target *KCNJ2*.

## References

 Krul SP, Driessen AH, van Boven WJ, et al. Thoracoscopic video-assisted pulmonary vein antrum isolation, ganglionated plexus ablation, and periprocedural confirmation of ablation lesions: first results of a hybrid surgicalelectrophysiological approach for atrial fibrillation. Circ Arrhythm Electrophysiol 2011;4:262–270.
Viviano A, Yin X, Zampetaki A, et al. Proteomics of the epicardial fat secretome and its role in post-operative atrial fibrillation. Europace 2018; 20:1201–1208.

3. Eyileten C, Jakubik D, Shahzadi A, et al. Diagnostic performance of circulating

miRNAs and extracellular vesicles in acute ischemic stroke. Int J Mol Sci 2022;23:4530.

4. Buntsma NC, Ga, secka A, Roos Y, et al. EDTA stabilizes the concentration of platelet-derived extracellular vesicles during blood collection and handling. Platelets 2022;33:764–771.

5. Welsh JA, Van Der Pol E, Arkesteijn GJA, et al. MIFlowCyt-EV: a framework for standardized reporting of extracellular vesicle flow cytometry experiments. J Extracell Vesicles 2020;9:1713526.

6. Gasecka A, Nieuwland R, Budnik M, et al. Ticagrelor attenuates the increase of extracellular vesicle concentrations in plasma after acute myocardial infarction compared with clopidogrel. J Thromb Haemost 2020;18:609–623.

7. Ruijter JM, Ramakers C, Hoogaars WM, et al. Amplification efficiency: linking baseline and bias in the analysis of quantitative PCR data. Nucleic Acids Res 2009;37:e45.

 PotseM, Linnenbank AC, Grimbergen CA. Software design for analysis of multichannel intracardial and body surface electrocardiograms. Comput Methods Programs Biomed 2002;69:225–236.

9. Lammers WJ, Schalij MJ, Kirchhof CJ, Allessie MA. Quantification of spatial inhomogeneity

in conduction and initiation of reentrant atrial arrhythmias. Am J

Physiol 1990;259:H1254-H1263.

10. van Eijndhoven MAJ, Aparicio-Puerta E, Gómez-Martín C, et al. Unbiased and UMI-informed sequencing of cell-free miRNAs at single-nucleotide resolution. bioRxiv https://doi.org/10.1101/2021.05.04.442244.

11. Aparicio-Puerta E, Gómez-Martín C, Giannoukakos S, et al. sRNAbench

and sRNAtoolbox 2022 update: accurate miRNA and sncRNA profiling for model and non-model organisms. Nucleic Acids Res 2022;50:

W710–W717.

12. Kozomara A, Birgaoanu M, Griffiths-Jones S. miRBase: from microRNA sequences to function. Nucleic Acids Res 2018;47:D155–D162.

13. Aparicio-Puerta E, Gómez-Martín C, Giannoukakos S, et al. mirnaQC: a webserver for comparative quality control of miRNA-seq data. Nucleic Acids Res


2020;48:W262–W267.

14. Gerstner N, Kehl T, Lenhof K, et al. GeneTrail 3: advanced high-throughput enrichment analysis. Nucleic Acids Res 2020;48:W515–W520.

15. McGeary SE, Lin KS, Shi CY, et al. The biochemical basis of microRNA targeting efficacy. Science 2019;366:eaav1741.


16. Chen Y, Wang X. miRDB: an online database for prediction of functional micro-RNA Targets. Nucleic Acids Res 2020;48:D127-D131.

## **Supplemental Figures and Tables**



### Supplemental Figure 1. Secretome samples concentrations.


Concentration of EAT and SAT secretome samples used for Flow cytometry or miRNA sequencing experiments expressed in mg of adipose tissue per mL of secretome. There was no significant differences in concentrations between the groups. One-way ANOVA. n=7 minimum.



Supplemental Figure 2. miR-133a-3p overexpression reduces *Kcnj12* expression.

Relative expression of *Kcnj12*, gene encoding for Kir2.2, in NRVMs transfected with NC and miR-133a-3p, measured by real-time qPCR. Data are mean±SD. NC: n=10, miR-133a-3p:

n=10 from three independent NRVMs isolation. Student's t-test. \**P*<.05. Primer sequences are included in Supplemental Table 4.



## Supplemental Figure 3. Ion channels expression

Relative expression of *Kcnj11*, *Kcnq1*, *Cacna1c*, *Cacna1g*, *Scn5a*, *Scn1b*, *Scn2b*, *Scn3b*, *Slc9a1*, *Atp1a1* and *Atp1a2* in NRVMs transfected with siNC, miR-1-3p and miR-133a-3p, measured by real-time qPCR. Data are mean±SD. siNC: n=10, miR-1-3p: n=10, miR-133a-3p: n=10 from three independent NRVMs isolation. One-way ANOVA followed by Tukey's multiple comparisons test in case of normal distribution, or by Dunn's multiple comparison test. \*P<.05. \*\*P<.01. Primer sequences are included in Supplemental Table 4.

**KCNJ2** Human 2729-GGCUUUUCCUU--UUGCACAUUCCA-2753

# *Kcnj2* Rat 2866-AGCUUU**U**CUU**U**CUU**U**GC**ACAUUCCA**-2890

KCNJ2 Human2729-GGCUUUUCCUU-UUGCACAUUCCA-2753hsa-miR-1-3p3'-UAUGUAUGAA-GAAUGUAAGGU-5'

<i>Kcnj2</i> Rat	2866–AGCUUU <mark>U</mark> CUU <mark>U</mark> CUU <mark>U</mark> GC <mark>ACAUUCCA</mark> –2890
hsa-miR-1-3p	3'-UAUGUAUGAAGAAAUGUAAGGU-5'

### Supplemental Figure 4. miR-1-3p predicted binding sites on KCNJ2

Complementary sequences between miR-1-3p and the predicted 3'UTRs binding sites of *KCNJ2* messenger RNA in human (ENST00000243457.4, ensembl.org) and rat (ENSRNOG00000064933, ensembl.org).

KCNJ12 Human	2430-AUUGCUUGUGCUAAA <mark>GACCAAA</mark> A-2454
<i>Kcnj12</i> Rat	48406-AUUGCUCGUGCUAAAGACCAAAC-48429
KCNJ12 Human	2430-AUUGCUUGUGCUAAA <mark>GACCAAA</mark> A-2454 
hsa-miR-133a-3p	3'- GUCGACCAACUUCCCCCUGGUUU -5'
Kcnj12 Rat	48406-AUUGCUCGUGCUAAA <mark>GACCAAA</mark> C-48429
hsa-miR-133a-3p	3'- GUCGACCAACUUCCCCUGGUUU -5'

### Supplemental Figure 5. miR-133a-3p predicted binding sites on *KCNJ12*

Complementary sequences between miR-1-3p and the predicted 3'UTRs binding sites of *KCNJ12* messenger RNA in human (ENST00000583088.6, ensembl.org) and rat (ENSRNOG0000002303, ensembl.org). Target sequence identified using targetscan.org.

# Supplemental Table 1. Patient characteristics

Experiment	All experiments	All experiments		y (paired secretome	miRNA sequencing		miRNA validation in EAT and SAT	
Clinical characteristics	n, mean±SD (range)	%	n, mean±SD (range)	%	n, mean±SD (range)	%	n, mean±SD (range)	%
Total	17		7		8		8	
Sex								
Male	10	58.8%	6	85.7%	5	62.5%	5	62.5%
Female	7	41.2%	1	14.3%	3	37.5%	3	37.5%
Age (years)	63±8 (39-75)		63±5 (57-73)		64±6 (57-75)		61±9 (39-70)	
Body Mass index (kg/m <sup>2</sup> )	28±3 (24-36)		28±2 (25-33)		27±3 (24-34)		29±4 (25-36)	
Atrial Fibrillation	× ,				· · ·			
Persistent AF	16	94%	7	100%	8	100%	7	87.5%
Long-standing persistent AF	1	6%	0	0%	0	0%	1	12.5%
Comorbidities								
Hypertension	10	59%	7	100%	6	75%	3	38%
Vascular diseases	0	0%	0	0%	0	0%	0	0%
Diabetes	2	12%	0	0%	2	25%	0	0%
Previous stroke	2	12%	1	14%	1	13%	0	0%
Heart failure	4	24%	1	14%	3	38%	1	13%
Age groups								
Age $\geq 75$	1	6%	0	0%	1	13%	0	0%
Age 65-74	6	35%	2	29%	2	25%	3	38%
Age <65	10	59%	5	71%	5	63%	5	63%

# Supplemental Table 2

miRNA	Averaged read-counts
hsa-let-7b-5p	231807
hsa-let-7a-5p	166820
hsa-let-7c-5p	138335
hsa-miR-26a-5p	52645
hsa-let-7f-5p	43885
hsa-miR-99a-5p	32951
hsa-miR-125b-5p	23020
hsa-miR-100-5p	17556
hsa-miR-320a-3p	16617
hsa-miR-143-3p	13779
hsa-miR-10b-5p	13681
hsa-miR-92a-3p	13441
hsa-let-7e-5p	12678
hsa-miR-30d-5p	11607
hsa-miR-148a-3p	9971
hsa-miR-125a-5p	8556
hsa-miR-1-3p	8378
hsa-miR-26b-5p	7208
hsa-miR-10a-5p	7020
hsa-miR-199a-3p	6766
hsa-miR-21-5p	6148
hsa-miR-23a-3p	5949
hsa-miR-101-3p	5640
hsa-miR-193a-5p	5586
hsa-miR-99b-5p	5543
hsa-miR-25-3p	5520
hsa-miR-486-5p	5268
hsa-miR-423-5p	5046
hsa-let-7i-5p	4379
hsa-miR-146b-5p	4239
hsa-miR-151a-3p	3970
hsa-miR-152-3p	3774
hsa-miR-574-5p	3770
hsa-miR-214-3p	3738
hsa-let-7d-5p	3700
hsa-miR-199b-3p	3347
hsa-miR-92b-3p	3300

hsa-miR-30a-5p	3079
hsa-miR-218-5p	3068
hsa-let-7g-5p	2893
hsa-miR-28-3p	2540
hsa-miR-222-3p	2462
hsa-miR-191-5p	2087
hsa-miR-127-3p	2060
hsa-miR-30e-3p	2035
hsa-miR-532-5p	1795
hsa-miR-3960	1463
hsa-miR-30e-5p	1453
hsa-miR-423-3p	1316
hsa-miR-342-3p	1300
hsa-miR-23b-3p	1249
hsa-miR-30a-3p	1196
hsa-miR-140-3p	1102
hsa-miR-24-3p	1095
hsa-let-7d-3p	1094
hsa-miR-451a	1092
hsa-miR-27a-3p	1079
hsa-miR-30c-5p	1078
hsa-miR-150-5p	1055
hsa-miR-103a-3p	1041
hsa-miR-574-3p	956
hsa-miR-27b-3p	949
hsa-miR-98-5p	915
hsa-miR-29a-3p	882
hsa-miR-328-3p	811
hsa-miR-432-5p	784
hsa-miR-181a-5p	760
hsa-miR-378a-3p	734
hsa-miR-1246	699
hsa-miR-374b-5p	693
hsa-miR-139-5p	665
hsa-miR-148b-3p	657
hsa-miR-126-3p	651
hsa-miR-184	633
hsa-miR-221-3p	629
hsa-miR-361-3p	628
hsa-miR-4510	616
hsa-miR-374a-5p	611
· · · · ·	

hsa-miR-195-5p	606
hsa-miR-199b-5p	576
hsa-miR-382-5p	566
hsa-miR-181b-5p	562
hsa-miR-1307-3p	559
hsa-miR-155-5p	544
hsa-miR-320b	520
hsa-miR-409-3p	507
hsa-miR-224-5p	500
hsa-miR-133a-3p	489
hsa-miR-22-3p	483
hsa-miR-126-5p	470
hsa-miR-125b-2-3p	460
hsa-miR-197-3p	442
hsa-miR-335-3p	432
hsa-miR-193b-5p	424
hsa-miR-379-5p	417
hsa-miR-4516	386
hsa-miR-186-5p	385
hsa-miR-4488	382
hsa-miR-93-5p	381
hsa-miR-361-5p	357
hsa-miR-485-5p	350
hsa-miR-204-5p	349
hsa-miR-146a-5p	340
hsa-miR-30b-5p	335
hsa-miR-483-5p	328
hsa-miR-132-3p	312
hsa-miR-151a-5p	307
hsa-miR-128-3p	302
hsa-miR-223-3p	290
hsa-miR-484	285
hsa-miR-199a-5p	283
hsa-miR-34a-5p	258
hsa-miR-1180-3p	240
hsa-miR-134-5p	238
hsa-miR-190a-5p	235
hsa-miR-320c	234
hsa-let-7b-3p	233
hsa-miR-744-5p	233
hsa-miR-106b-3p	226
<b>•</b>	

hsa-miR-1260a215hsa-miR-99b-3p204hsa-miR-4454204hsa-miR-16-5p200hsa-miR-874-3p195hsa-miR-874-3p195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-110175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p143hsa-miR-660-5p143hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p106hsa-miR-877-5p101	hsa-miR-203a-3p	217
hsa-miR-99b-3p204hsa-miR-4454204hsa-miR-16-5p200hsa-miR-874-3p195hsa-miR-320d195hsa-miR-320d195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p143hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-671-5p127hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p109		
hsa-miR-4454204hsa-miR-16-5p200hsa-miR-874-3p195hsa-miR-320d195hsa-miR-320d195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-1247-5p175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-340-5p109		
hsa-miR-16-5p200hsa-miR-874-3p195hsa-miR-320d195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-1247-5p175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-135a-5p157hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-660-5p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p106hsa-miR-340-5p101		
hsa-miR-874-3p195hsa-miR-320d195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-92b-5p157hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p127hsa-miR-660-5p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p109hsa-miR-340-5p109hsa-miR-340-5p109		
hsa-miR-320d195hsa-miR-1247-5p184hsa-miR-1247-5p184hsa-miR-2110175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-660-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p109hsa-miR-340-5p109hsa-miR-340-5p109		
hsa-miR-1247-5p184hsa-miR-2110175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-660-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p101		
hsa-miR-2110175hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-660-5p127hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p109		
hsa-miR-139-3p173hsa-miR-135a-5p171hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-660-5p143hsa-miR-660-5p143hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365b-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p101	<b>L</b>	
hsa-miR-135a-5p171hsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p143hsa-miR-660-5p143hsa-miR-660-5p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-365a-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p101		
Image: Non-relationImage: Non-relationhsa-miR-543165hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-342-5p154hsa-miR-704150hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-671-5p127hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-92b-5p157hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-342-5p154hsa-miR-7704150hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-660-5p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101	-	
hsa-miR-1260b155hsa-miR-342-5p154hsa-miR-340-3p144hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-365a-3p109hsa-miR-340-5p109hsa-miR-340-5p109hsa-miR-340-5p101		
hsa-miR-342-5p154hsa-miR-7704150hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-7704150hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-340-3p144hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-505-3p127hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101	-	
hsa-miR-452-5p143hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-660-5p143hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-505-3p138hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-671-5p127hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101	-	
hsa-miR-365b-3p121hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-365a-3p121hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-369-5p118hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		
hsa-miR-1301-3p109hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101		121
hsa-miR-340-5p109hsa-miR-433-3p106hsa-miR-877-5p101	hsa-miR-369-5p	118
hsa-miR-433-3p106hsa-miR-877-5p101	hsa-miR-1301-3p	109
hsa-miR-877-5p 101	hsa-miR-340-5p	109
-	hsa-miR-433-3p	106
has miD 1914 5m 100	hsa-miR-877-5p	101
nsa-mik-1810-5p 100	hsa-miR-181d-5p	100
hsa-miR-365a-5p 100	hsa-miR-365a-5p	100
hsa-miR-185-5p 99	hsa-miR-185-5p	99
hsa-miR-193b-3p 90	hsa-miR-193b-3p	90
hsa-miR-425-5p 80	hsa-miR-425-5p	80
hsa-miR-374a-3p 67	hsa-miR-374a-3p	67
hsa-miR-145-3p 66	hsa-miR-145-3p	66
hsa-miR-1290 63	hsa-miR-1290	63
hsa-miR-1306-5p 53	hsa-miR-1306-5p	53

# Supplemental Table 3. miRNA mimic sequences

miRNA mimic	Sequence 5'-3'
hsa-miR-1-3p	UGGAAUGUAAAGAAGUAUGUAU
hsa-miR-133a-3p	UUUGGUCCCCUUCAACCAGCUG

# Supplemental Table 4. Primer sequences

Primer	Sequence					
Atp1a1 Fw	CGGATGACAGTGGCTCACAT					
Atp1a1 Rv	CGCTACTGCACGCTTAAGGA					
Atp1a2 Fw	CATCATTTGCAAGACCCGGC					
Atp1a2 Rv	AACCACCACGTGACCTTGAG					
Cacna1c Fw	TGGGATCATGGCTTATGGCGGC					
Cacna1c Rv	ATCAGCCAGGTTGTCCACCG					
Cacna1g Fw	AGGCAGAGGAAATCGGCAAA					
Cacna1g Rv	CTGTCCCCATCACCATCCAC					
Eef1e1 Fw	TCCAGTAAAGAAGACACCCAGA					
Eef1e1 Rv	GACAAAACCAGCGAGACACA					
Gapdh Fw	GGTGGACCTCATGGCCTACA					
Gapdh Rv	CTCTCTTGCTCTCAGTATCCTTGCT					
Hprt Fw	TGACTATAATGAGCACTTCAGGGATTT					
Hprt Rv	CGCTGTCTTTTAGGCTTTGTACTTG					
Kcnj2 Fw	TGTGTTACAGACGAGTGCCC					
Kcnj2 Rv	CAGAGTTTGCCGTCCCTCAT					
Kcnj11 Fw	ATCAGTCCAGAGGTTGGTGC					

V 111 D	
Kcnj11 Rv	TAATGCCCTTTCGGGACAGC
Kcnj12 Fw	AACCCCTACAGCATCGTATC
Kcnj12 Rv	GCACCTTGCCATTGCCAAA
Kcnq1 Fw	GATCAGTCCATCGGGAAGCC
Kcnq1 Rv	GGTCCAGTTGTGTCACCTTGT
Slc9a1 Fw	AACGGCTGCGGTCCTATAAC
Slc9a1 Rv	CGAGACATGGTGGGTGAGTC
Scn1b Fw	AACACCAGCGTCGTCAAGAA
Scn1b Rv	TTCCGAGGCATTCTCTTGTGC
Scn2b Fw	CCTTGGTCCCTCAATCACCC
Scn2b Rv	ACTGTGACTTCCATGCTCCG
Scn3b Fw	CTTCCTCACCTTGTGGCTGT
Scn3b Rv	GCACTCAGATCACCTCAAGTCA
Scn5a Fw	TCTTCCGGTTCAGTGCCACC
Scn5a Rv	GGATGGTGCACATGATGAGCATG
U6 Fw	CTCGCTTCGGCACA
U6 Rv	AACGCTTCACGAATTTGCGT
hsa-miR-1-3p	TGGAATGTAAAGAAGTATGTAT
hsa-miR-133a-3p	TTTGGTCCCCTTCAACCAGCTG
hsa-miR-195-5p	TAGCAGCACAGAAATATTGGC
hsa-miR-24-3p	TGGCTCAGTTCAGCAGGAACAG
hsa-miR-26a-5p	TTCAAGTAATCCAGGATAGGCT
hsa-miR-26b-5p	TTCAAGTAATTCAGGATAGGT
L	1

hsa-miR-328-3p	CTGGCCCTCTCTGCCCTTCCGT

# MIFlowCyt-EV of "MicroRNAs in Extracellular Vesicles released from Epicardial Adipose Tissue Promote Arrhythmogenic Conduction Slowing"

# **Flow cytometry**

# **Experimental design**

The aim of flow cytometry (A60-Micro, Apogee Flow Systems, Hemel Hempstead, UK) experiments was to measure the concentrations of extracellular vesicles (EVs) released from epicardial and subcutaneous adipose tissue (EAT, SAT) in conditioned culture medium from adipose tissue samples collected from 7 patients to investigate whether EAT and SAT secretome (EATs, SATs) show different concentrations of EVs. We hypothesized that EAT secretes a larger number of EVs than SAT.

All samples were measured using an autosampler, which facilitates subsequent measurements of samples in a 96-well plate. For this study, samples were measured on two different days (03/01/2021 for 3 patient samples and 04/05/2022 for 4 patient samples). Each day a buffer-only control, antibody in buffer controls, and isotype controls corresponding to the labels in the well plate were measured. Flow rate and scatter calibrations were performed daily. Fluorescence calibration was performed twice during this time period. To automatically determine optimal sample dilutions, apply calibrations, determine and apply gates, generate reports with scatter plots and generate data summaries, we used custom-built software (MATLAB R2020b, Mathworks, Natick, MA, USA).

## Sample dilutions

Samples containing submicrometer particles require different dilutions to (1) avoid swarm detection [1] and (2) achieve statistically significant counts within a clinically applicable measurement time. Although serial dilutions are recommended to find the optimal dilution, we consider serial dilutions unfeasible in clinical research studies. Therefore, we developed a procedure to estimate to optimal sample dilution (section 1.2 of <u>https://doi.org/10.6084/m9.figshare.c.4753676</u>). In sum, we showed that for our flow cytometer and settings used, a count rate  $\leq 5.0 \cdot 10^3$  events second unlikely results in swarm detection.

Pre-staining control and SATs samples were undiluted, while EATs samples were undiluted, 1.5x diluted or 2x diluted, followed by an additional 11.125x dilution post-

staining, resulting in count rates between 18 and 3844 events per second. More research is required, however, to confirm that the count rate can be used as a benchmark to avoid swarm detection for a given sample type.

# **EV staining**

EVs were stained with antibodies. Prior to staining, antibodies were diluted in DPBS and centrifuged at 18,890 g for 5 min to remove aggregates. **Supplemental Table 5** shows an overview of the used reagents and antibody concentrations during staining. Each sample was single labelled with CD9-PE, CD63-PE or Lactadherin-FITC. To stain, 20  $\mu$ L of pre-staining diluted EATs and SATs samples was incubated with 2.5  $\mu$ L of antibodies or isotype controls and kept in the dark for 2 h at room temperature. Post-staining, samples were diluted 11.125-fold in 200  $\mu$ L of DPBS to decrease background fluorescence from unbound reagents.

## **Buffer-only control**

Each measurement day at least 1 well with DPBS was measured with the same flow cytometer and acquisition settings as the samples. The mean count rate was 16 events per second, which is substantially lower than the count rates obtained for EATs/SATs samples (18 - 3,844 events per second).

## **Buffer with reagents control**

Each 96-well-plate contained a buffer with reagent control for each reagent (**Supplemental Table 5**), which was measured with the same flow cytometer and acquisition settings as the samples. For all reagents, the mean count rate was between 14 - 22 events per second, which is in the same range as the buffer-only control.

# **Unstained controls**

Unstained controls were measured with the same flow cytometer and acquisition settings as the stained samples, resulting in 15 - 3,427 events per second for EATs/SATs samples.

## **Isotype controls**

**Supplemental Table 5** shows an overview of the used isotype controls, which were added to a selection of samples. For 21 control/EATs/SATs samples, we obtained an average of  $5.5 \times 10^4$  IgG1-PE+ events with a diameter 160 nm - 1,000 nm per measurement. For comparison, the average concentration was  $2.61 \times 10^5$  CD63-PE+,  $4.46 \times 10^5$  CD9-PE+,  $5.20 \times 10^6$  Lactadherin-FITC+ events per mL with a diameter 160 nm - 1,000 nm were obtained in the experiments using EATs/SATs samples.

# **Detergent-lysis controls**

20 μL of 10% NP-40 Detergent Solution were added to 2 of the EATs and 2 of the SATs labelled sample for detergents control analysis. Results are provided in **Supplemental Figure 6** and show that the detergent step removes EV signal in EATs and SATs.

# Trigger channel and threshold

Based on the buffer-only control (38 events s<sup>-1</sup>), the acquisition software was set up to trigger at 14 or 24 arbitrary units (a.u.) SSC, which is equivalent to an SSC cross section of 8 or 10 nm<sup>2</sup> (Rosetta Calibration, v1.29, Exometry, Amsterdam, The Netherlands), respectively. Both trigger thresholds enable detection of EVs with a diameter of 160 nm, which is the lower boundary of one of the applied gates.

# Flow rate quantification

The A60-Micro is equipped with a syringe pump with volumetric control, which we checked on a daily basis using Apogee Mix and ApoCal Mix. Here, we assumed that the flow rate is equal to the adjusted flow rate of  $3.01 \,\mu$ L/min for all measurements.

# **Fluorescence** calibration

Calibration of the fluorescence detectors from arbitrary units (a.u.) to molecules of equivalent soluble fluorochrome (MESF) was accomplished using SPHERO Easy Calibration Fluorescent Particles (AK01, Spherotech Inc., Irma Lee Circle, IL, USA) and Quantum<sup>TM</sup> MESF Kits (13734, Bangs Laboratories Inc., Fishers, IN, USA).

Calibration of the FITC and PE detector was performed once in the time period of these experiments. For each measurement, we added fluorescent intensities in MESF to the flow cytometry data files by custom-built software (MATLAB R2020b) using the following equation:

$$I(MESF) = 10^{a \cdot \log_{10} I(a.u.) + b}$$
Equation S1

where I is the fluorescence intensity, and *a* and *b* are the slope and the intercept of the linear fits, respectively, see **Supplemental Table 6**.

# Light scatter calibration

We used Rosetta Calibration to relate scatter measured by FSC or SSC to the scattering cross section and diameter of EVs. Fig. S7 shows print screens of the scatter calibrations. We

modelled EVs as core-shell particles with a core refractive index of 1.38, a shell refractive index of 1.48, and a shell thickness of 6 nm.

# MIFlowCyt checklist

The MIFlowCyt checklist is added to **Supplemental Table 7**.


# EV number concentration

The concentrations reported in the manuscript describe the number of particles (1) that exceeded the SSC threshold, corresponding to a side scattering cross section of 8 nm<sup>2</sup> or 10 nm<sup>2</sup>, (2) that were collected during time intervals, for which the count rate did not a deviated more than 750 events/s of the median count rate, (3) with a diameter 160 nm – 1,000 nm as measured by SSC after light scatter calibration (section 1.11) and (4) are positive for PE or FITC, per mL of sample. The lower fluorescence gate is 148 MESF for CD63-PE, 130 MESF for CD9-PE, and 512 MESF for Lactadherin-FITC.

# Data sharing


Data is available via:

 $https://figshare.com/articles/dataset/Flow\_cytometry\_of\_Epicardial\_adipose\_tissue\_derived-extracellular\_vesicles\_/21695444$ 



Supplemental Figure 6. Detergent step removes EV signal in EATs and SATs.

Concentrations of EVs per mL in EATs (A) and SATs (B) from n=2 patients, measured by flow cytometry after labelling EVs with anti-CD9-PE, anti-CD63-PE, and anti-Lactadherin-FITC, with and without detergent step (cf methods). Concentrations of EVs show the number of particles per mL of secretome (i) exceeding the side scattering threshold (ii) having a diameter <1,000 nm, and (iii) exceeding the fluorescent gate corresponding to the used labels. Data are mean  $\pm$  standard error of the mean.



**Supplemental Figure 7. Scatter calibrations.** (A) Forward scatter and (B) side scatter calibration of the A60-Micro by Rosetta Calibration. To relate scatter to the diameter of EVs, we modelled EVs as core-shell particles with a core refractive index of 1.38, a shell refractive index of 1.48, and a shell thickness of 6 nm.

**Supplemental Table 5: Overview of staining reagents.** Characteristics being measured, analyte, analyte detector, reporter, isotype, clone, concentration, manufacturer, catalog number and lot number of used staining reagents. The antibody concentration during measurements was 11.25-fold lower than the antibody concentration during staining.

Characteristic measured	Analyte	Analyte detector	Reporter	Isotype	Clone	Concentration during staining (µg mL <sup>-1</sup> )	Manufacturer	Catalog number	Lot number
Tetraspanin	Human CD9	Anti-human CD9 antibody	PE	IgG1	M-L13	2.5	BD bioscience	555372	03905
Tetraspanin	Human CD63	Anti-human CD62 antibody	PE	IgG1	CLBGran/ 12	7.5	Beckman Coulter	IM1914	28
Phosphatidyl-L- serine	Bovine Lactadherin	Phosphatidyl-L-serine	FITC	IgG1	n.a.	10.4	Prolytix	BLAC- FITC	MM0301
Affinity for Fc receptor	Fc receptor	IgG1	PE	n.a.	X40	7.5	BD Bioscience	345816	9309643

CD: cluster of differentiation; FITC: fluorescein isothiocyanate; IgG: immunoglobulin G; PE: phycoerythrin.

	Calibration date	Slope	Intercept	R <sup>2</sup>
FITC	2020-07-20	1.2291	-1.4831	0.9988
FITC	2022-03-30	1.1823	-2.1922	0.9993
PE	2020-07-20	1.0340	-1.6142	0.9993
PE	2022-03-29	1.2538	-2.7574	0.9916

# Supplemental Table 6: Overview of fluorescence calibrations.

FITC: fluorescein isothiocyanate; PE: phycoerythrin. R<sup>2</sup>: coefficient of determination.

# Supplemental Table 7 . MIFlowCyt checklist.

Requirement	Please Include Requested Information
1.1. Purpose	To compare the EVs content between EATs and SATs obtained
	from 7 patients with atrial fibrillation
1.2. Keywords	Epicardial adipose tissue, atrial fibrillation, extracellular vesicles,
	arrhythmias
1.3. Experiment variables	EAT and SAT paired samples were collected from 7 different
	patients, on different dates.
1.4. Organization name and	Amsterdam University Medical Centers
address	Location Academic Medical Centre
	Meibergdreef 9
	1105 AZ Amsterdam
	The Netherlands
1.5. Primary contact name	Auriane C. Ernault,
and email address	Amsterdam University Medical Centers
	Location Academic Medical Centre
	Meibergdreef 9
	1105 AZ Amsterdam
	The Netherlands
	a.c.ernault@amsterdamumc.nl
1.6. Date or time period of	Measurement was performed on 03/01/2021 for 3 patients
experiment	samples and 04/05/2022 for 4 patients samples
1.7. Conclusions	EATs contains significantly more EVs than SATs
1.8. Quality control measures	All samples were measured using an autosampler, which
	facilitates subsequent measurements of samples in a 96-well plate.
	Each well plate contained buffer-only controls (section S1.4),
	antibody in buffer controls (section S1.5), unstained controls
	(section S1.6) and isotype controls (section S1.7). The flow rate is
	calibrated with Apogee Calibration beads (Apogee Flow Systems,
	Hemel Hempstead, UK; section S1.9). Fluorescence detectors
	were calibrated (section S1.10) with SPHERO Easy Calibration
	Fluorescent Particles (AK01, Spherotech Inc., Irma Lee Circle, IL,
	USA) and Quantum <sup>TM</sup> MESF Kits (13734, Bangs Laboratories

Inc., Fishers, IN, USA). FSC and SSC were calibrated with
Rosetta Calibration (v1.29, section S1.11).
Samples were measured on 2 different days, in a time period of 17
months.
The samples analysed in this study are EATs and SATs.
Adipose tissue was rinsed in phosphate-buffered saline, dissected
to remove any blood or connective tissue, and minced into 1-
2mm <sup>3</sup> explants. Each explant piece was transferred to a well of a
96-well plate and incubated in 100 $\mu$ L of culture medium for 24H
at 37°C and 5% CO <sub>2</sub> in order to stimulate adipokine and
secretome release. Adipose tissue secretome was then collected.
For freeze-storage, secretome samples were transferred to 1.5 mL
cryotubes and stored in -80°C. Before staining, samples were
thawed for 1 minute at 37°C.
Human adipose tissue was collected as follow: seven patients with
symptomatic, persistent or long-standing persistent atrial
fibrillation (AF) undergoing video-assisted thoracoscopic
pulmonary vein isolation (VATS-PVI) were recruited at the
Amsterdam UMC, University of Amsterdam, between 2019 and
2022. Written informed consent was obtained from all
participants. The study was conducted in accordance with the
Declaration of Helsinki as revised in 2013 and with the Medical
Research Involving Human Subjects Act (WMO) and other local
laws and regulations. AF was documented on ECG, Holter or
pacemaker electrogram at least once in the 12 months preceding
surgery. Patients with left ventricular ejection fraction <30% and
severe heart failure were excluded from the study. During VATS-
PVI the left atrial appendage (LAA) was resected as standard
procedure. The atrial tissue was collected and the EAT overlying
the LAA was removed. A biopsy of approximately 1cm <sup>3</sup> of
subcutaneous adipose tissue (SAT) was obtained from one of the
thoracoscopic entry sites. Adipose tissue was transported in
phosphate-buffered saline to the laboratory.

2.2 Sample characteristics	Epicardial and subcutaneous adipose tissue secretome are
	expected to contain EVs, lipoproteins and proteins.
2.3. Sample treatment	Samples were thawed and brought to the lab by manual transport
description	in vertical position. Samples were prepared by double
	centrifugation and samples were stained. Please see section S1.3
	for staining procedure.
2.4. Fluorescence reagent(s)	Please see Supplemental Table 5.
description	
3.1. Instrument manufacturer	Apogee, Hemel Hempstead, UK
3.2. Instrument model	A60-Micro
3.3. Instrument configuration	Samples were analysed for 2 minutes at a flow rate of
and settings	3.01 $\mu$ L/min on an A60-Micro, equipped with a 405 nm laser
	(100 mW), 488 nm laser (150 mW) and 638 nm laser (150 mW).
	The trigger threshold was set at SSC 14 or 24 a.u. corresponding
	to an SSC cross section of 8 nm <sup>2</sup> or 10 nm <sup>2</sup> (Rosetta Calibration),
	respectively. For FSC, the PMT voltages were 380 V or 470 V.
	For SSC, the PMT voltages were 350 V or 375 V. For all
	detectors, the peak height was analysed. FITC signals were
	collected with the 488-Green(Peak) detector (525/50 nm band
	pass filter, PMT voltage 520 V and 560 V). PE signals were
	collected with the 488-Orange(Peak) detector (575/30 nm band
	pass filter, PMT voltage 520 V and 450 V).
4.1. List-mode data files	Data is available via FigShare.com: 10.6084/m9.figshare.21695444
4.2. Compensation	No compensation was required because no fluorophore
description	combinations were used that have overlapping emission spectra.
4.3. Data transformation	Fluorescence detectors were calibrated (section S1.10) with
details	SPHERO Easy Calibration Fluorescent Particles (AK01,
	Spherotech Inc., Irma Lee Circle, IL, USA) and $Quantum^{TM}$
	MESF Kits (15036, Bangs Laboratories Inc., Fishers, IN, USA).
	FSC and SSC were calibrated with Rosetta Calibration (v1.29,
	section S1.11). The concentrations reported in the manuscript

	describe the number of particles that fulfil the gating criteria per
	mL.
4.4.1. Gate description	To automatically apply gates, generate pdf reports with scatter
	plots, and summarize the data in a table, custom-built software
	(MATLAB R2018b) was used. Please find below a description of
	the gates. First, events that were collected during seconds for
	which the count rate deviated less than 750 events/seconds from
	the median count rate were included. Second, events with a
	diameter 160 nm - 1,000 nm as measured by SSC after light
	scatter calibration (section S1.11) were included. Third, events
	positive for either FITC or PE were included.
4.4.2. Gate statistics	The number of positive events was corrected for flow rate,
	measurement time and dilutions performed during sample
	preparation.
4.4.3. Gate boundaries	The lower boundary of the diameter gate was 160 nm. The upper
	boundary of the diameter gate was 1,000 nm. The lower
	boundaries of the fluorescent gates were automatically determined
	(MATLAB R2020b). The lower fluorescence gate is 148 MESF
	for CD63-PE, 130 MESF for CD9-PE, and 512 MESF for
	Lactadherin-FITC.

a.u.: arbitrary units; EVs: extracellular vesicles; FSC: forward scattering; SSC: side scattering.